

# Command Set CenSon

0,6

## COM Settings:

Baudrate: 9600  
 Databit: 8  
 Parity: none  
 Stop Bit 1  
 Handshake: OFF

From PC to Device	Answer
-------------------	--------

CMD	Frame	Params	Unit	Range	Ack1	Ack2
<b>General</b>						
#DEVINI_T0	<CR>	-	-	-	Ack-	INI
<b>Sonicator (positioning):</b>						
#SONPOS_W0	<CR>	-	-	-	Ack-	SPW
#SONPOS_Sn	<CR>	n=[1..100]	mm	1-100mm	Ack-	SPS
<b>Sonicator (sonification):</b>						
#SONSNC_Pn	<CR>	n=[1..255]	0.1 sec	0.1-25.5s	Ack-	SSP
#SONSNC_Wn	<CR>	n=[1..255]	0.1 sec	0.1-25.5s	Ack-	SSW
<b>Centrifuge (positioning):</b>						
#CENPOS_P1	<CR>	-	-	-	Ack-	CP1
#CENPOS_P2	<CR>	-	-	-	Ack-	CP2
<b>Centrifuge (centrifugation):</b>						
#CENSET_Sn	<CR>	n=[10..100]	10 rot/min	100-1000rot	Ack-	CSS
#CENRUN_Tn	<CR>	n=[1..180]	10 sec	10s - 30min	Ack-	CRUN
#CENSTA_T0	<CR>	-	-	-	Ack-	CRDY
						BUSY
#CENRUN_Dn	<CR>	n=[1..180]	10 sec	10s - 30min	Ack-	
else					Ack-	Err

Info
------

Timeout	Action	Behaviour
60s	Washer: ~6sec, Centrifuge: ~8sec	
30s	Move arm to wash pos.	
30s	Move arm to sonic pos.first, then dive-in n [mm]	locked
30s	Sonicate in probe for [n] 1/10seconds	unlocked
30s	SONSNC_W10 = WashPump 1sec	unlocked
20s	Present position 1	locked
20s	Present position 2	locked
10s	Set rotation speed in steps of 100 [100, 200, ... 1000]	
300s	Start centrifugation for [n] seconds	
no	Centrifuge is ready for tube transfers (safe)	
	Centrifuge is busy with CENRUN_Dn	
no	No timeout, no second acknowledge	unlocked
	Any mismatch in commands returns Err	

## Start Behaviour

- After power on, the device does automatically init the sonic arm and moves the centrifuge to position 1
- The Sonic arm should always be centered over the wash station before power on the system! (ZeroPoint reference)
- Starting and initializing EVOware sends automatically a init command to all devices, in this case the command: DEVINI\_T0

## Timeout

- Timeouts are handled on the PIPMOD driver side. If no matching acknowledge combination can be caught within the specified timeout, then the driver will report an error to EVOware. Not yet for all commands

## Command Behaviour

- After sending a command, there is a immediate response with Ack1 (see list above)
- When the command is done, a second Acknowledge is sent: Ack2
- The controller is quiet during this time, no statuscheck or other commands can be used until the second acknowledge is finished
- EXCEPTION: the SONSNC\_xx commands are unlocked, so the z pos can be changed during sonification

## Revision 0.6

- new command added to allow centrifuge to run parallel
- Status command extended med BUSY